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
Stereotactic Radiosurgery Comparison			
Features	Cyberknife	GammaKnife	TomoTherapy
Number of beams	50 - 300	≤ 201	1,000 - 50,000
Angles	Non-coplanar	Non-coplanar	Coplanar
Total treatment time	30 - 180 min	40 - 70 min	30 - 75 min
Field size	5 - 60mm diameter	4 - 18mm diameter	6mm x 1 > 6mm x 5
Delivery	Fixed, circular	Fixed, circular	Helical
Positioning accuracy	sub - mm	sub - mm	sub - mm
Energy	6MV	~1.25MV, Cobalt	6MV
Tx room construction	Large footprint	Small	Small (20 x 15 x 9 ft)
Shielding	Primary & Secondary	Secondary	Secondary
Clinical range	Cranial/Extracranial	Cranial	Cranial/Extracranial
	SRS / SBRT	SRS	SRS / SBRT & IMRT
Reimbursement	1 - 5 fractions	1 fraction	1 - 5, or IMRT > 40
Immobilization	Frameless	Frame	Frame or Frameless
Tracking	Yes	No	No
Positioning	2D x-ray	External Coord.	3D CT
Delivery verification	No	No	Dose calc, 3DCT
Staffing	Dedicated	Dedicated	Existing RT
Installation to clinical	TBD	?	< 30 days

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
Soisson ET, Hoban PW, Kammeyer T, Kapatoes JM, Westerly DC, Basavatia A, Tomé WA. A technique for stereotactic radiosurgery treatment planning with helical tomotherapy. *Med Dosim.* 2010 Feb 4.

Soisson, E.T.; Sobering, G.; Lucas, D.; et al. Quality assurance of an image guided intracranial stereotactic positioning system. *Technol. Cancer Res. Treat.* 8:39-49; 2009.

Bichay, T.; Chen, C.; Davis, S.; et al. TomoTherapy for cranial radiosurgery/radiotherapy. *IFMBE Proc.* 14:1852-5; 2007.

Holmes, T.W.; Hudes, R.; Dziuba, S.; et al. Stereotactic image guided intensity modulated radiotherapy using the Hi-ART II helical tomotherapy system. *Med. Dosim.* 33:135-48; 2008.

Yartsev, Y.; Kron, T.; Cozzi, L.; et al. Tomotherapy planning of small brain tumors. *Radiother. Oncol.* 74:49-52; 2005.

Han, C.; Liu, A.; Schultheiss, T.W.; et al. Dosimetric comparisons of helical tomotherapy treatment plans and step-and-shoot intensity-modulated radiosurgery treatment plans in intracranial stereotactic radiosurgery. *Int. J. Radiat. Oncol. Biol. Phys.* 65:608-16; 2006.

Gutiérrez, A.; Tomé, WA.; Stathakis, S.; et al. Integrated subvolume boosting technique to reduce normal Tissue irradiation during SBRT treatment planning using helical tomotherapy. *Med. Phys.* 25:2957; 2008.

TomoTherapy Incorporated. Field Safety Notice 3449. 2008. Available at: http://tomoexchange.tomotherapy.com/documentation/category/product_advisory_notices. Accessed December 16, 2009.

Woch, K.; Soisson, E.; Westerly, D.; et al. Acceleration of tomotherapy treatment delivery by increasing pitch and decreasing modulation. *Med. Phys.* 35:2636-7; 2008.

Kissick, M.W.; Flynn, R.T.; Westerly, D.C.; et al. On the making of sharp longitudinal dose profiles with helical tomotherapy. *Phys. Med. Biol.* 52:6497-510; 2007.